

yuwell



yuwell

JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Manufacturing Address: No.1 Baisheng Road, Development Zone, Danyang, Jiangsu Province, P.R.China
Tel: 86-0511-8690 0859 86-0511-8690 0860
E-mail: yuwellsales@yuyue.com.cn; sales@yuyue.com.cn
www.yuwell.com

f | X | in | @ | www.yuwell.com

2024, Yuyue medical, all rights reserved. The " yuyue ", yuwell, the yuwell logo and name are trademarked by Jiangsu Yuyue Medical Equipment & Supply Co., Ltd, and may not be distributed, modified, reproduced in whole or in part without the prior. All specifications and product configurations are subject to change without notification.

MEDICAL DEVICE
PRODUCT BROCHURE

CONTENTS

ABOUT YUWELL	01
RESPIRATORY	03
BreathCare PAP Series	03
Breathwear Mask Series	17
Oxygen Concentrators	21
Nebulizers	31
DIABETES MANAGEMENT	41
Blood Glucose Meters	41
Continuous Glucose Monitoring Systems	51
DIAGNOSTICS	55
Electronic Blood Pressure Monitors	55
Thermometers	71
Fingertip Pulse Oximeters	83
SUCTION	89
Portable Phlegm Suction Units	89
Suction Apparatus	95
REHABILITATION	101
Electric Wheelchairs	101
Mattresses with Mechanical Pump	109
YUWELL IOT PLATFORM	115

GLOBAL MANUFACTURING BASE

CHINA'S LARGEST HOME CARE
MEDICAL EQUIPMENT COMPANY

FOCUSED ON MEDICAL TREATMENT
FOR **26** YEARS


CPS-based smart devices
Artificial intelligence equipment leveraging visual recognition and detection
Intelligent warehousing and logistics system
Traceability platform leveraging big data, covering the entire life cycle of our medical devices

yuwell 鱼跃

GLOBAL DISTRIBUTION

12 R&D Centers / **9** Production Centers / **103** Countries & Regions

- Netherland
- Rottweil, Germany
- Tuttlingen, Germany
- Italy
- Beijing
- Shanghai
- Nanjing
- Suzhou
- Tibet
- Danyang
- Beijing
- Shanghai
- Danyang
- Suzhou
- Tibet
- Italy
- Germany



yuwell

BREATHCARE PAP SERIES

RESPIRATORY



NEW



Professional Care

Intelligent algorithm

The optimized therapy algorithm with additional respiratory-events monitoring guarantees a more comprehensive and comfortable treatment experience.

Ultra-quiet

The newly developed blower and airway design enable ultra-quiet operation of the BreathCare III.

Mask-fit test

A simple mask fitting test before sleep ensures the optimal fit, providing patients with greater confidence as they fall asleep.

Ergonomic Interaction

Convenient touch screen

The complex operations of the CPAP device have been simplified into easy-to-use touch controls, the large and responsive screen makes settings easier than ever.

Auto-all


Automized settings including auto-ramp, auto-humidity, auto-heated tube, auto-light, smart start/stop, altitude compensation and leak compensation, allowing the patients to receive treatment with minimum confusion and improved compliance.

Wireless data system

A seamless cloud-based system enables all patient care network participants to find their most suitable portal to therapy data, with convenient functionality.

Model	YH-680	YH-690
Picture		
Color	White/Black	
Certification	CE	CE
Mode	CPAP, Auto CPAP	CPAP, Auto CPAP
Pressure	4-20cmH ₂ O	4-20cmH ₂ O
Noise	27dBA±2dBA	27dBA±2dBA
Screen	3.5-inch	4.3-inch
Efficient Humidifier	•	•
FPS-Tech	•	•
Pre-heat		•
Auto-all	•	•
Central Apnea Detection	•	•
Leak Reminder	•	•
Data Management	•	•
Wireless Transfer	•	•
Cloud Platform	•	•
G.W./N.W.	3.7/2.4kg	
Packing	1SET/CTN 41×31×18cm	

INFORMATION INTERCONNECTION



INTELLIGENT DATA MANAGEMENT

Managing the therapy data through multiple platforms, seamlessly monitor the treatment process on your mobile phone, office computer or simply a web browser. Meeting various demands of different roles in the patient care network.



WIRELESS TRANFER

Comprehensive and professional respiratory therapy analysis system based on security-assured cloud server for patient management. Optional WiFi module or 4G Cellular module are required for this function. Currently available on BreathCare II, BreathCare Eco and BreathCare III

QR CODE

User-friendly sleep quality assessment with essential parameters. Unique QR code will be generated after each therapy and can be easily interpreted by scanning with Yuwell mobile app BreathCare+. Currently available on YH-360/YH-560 and YH-450/YH-480.


SD CARD/USB

Traditional and reliable data transfer method. Stored data can be interpreted with Yuwell computer software BreathCare Station and automatically generate comprehensive report. Available on all type of products.

YUWELL BREATHCARE WEB SYSTEM

It is a seamless, securely-hosted and cloud-based system for managing patients with sleep-disordered breathing. You can monitor your patients' device remotely without leaving your desk. The system makes it easier to simplify the workflow and to collaborate more efficiently across the patient care network.

Yuwell BreathCare Web System can now be associated with BreathCare II, BreathCare Eco and BreathCare III.



BREATHCARE ECO CPAP / AUTO CPAP



Stable and Reliable Quality

- High Performance Motor
Silent DC brushless motor, 0.2 seconds quick response, and 20000 hours of service life.
- No-Spoiler Air Duct Design
Accurately monitors the air flow signals and reduces noise at the same time.
- Patented Micro Airflow Laminar Flow Tube
Equipped with the world's advanced pressure and flow sensors. Accurately collects pressure & flow signal changes. Real-time control.

User-Friendly Design

- Convenient Button Operation
A simple combination of two buttons for all operations.
- Clear Treatment Information
One step to get all therapy information clearly on a 2.7 inch screen.
- Brand New Design Humidification Tank
Optimized humidification tank designed to better prevent water backflow.
- Yuwell Patient Management System
It is a seamless, securely-hosted and cloud-based system for managing patients with sleep-disordered breathing.

YUWELL BREATHCARE PAP / PRODUCT PARAMETERS

Specification	BreathCare Eco	
Model	YH-350	YH-550
Picture		
Color	White	
Certification	CE	CE
Mode	CPAP	CPAP, Auto CPAP
Pressure	4-20cmH ₂ O	4-20cmH ₂ O
Noise	28dBA±2dBA	28dBA±2dBA
Efficient Humidifier	•	•
Smart Start/Stop	•	•
Central Apnea Detection	•	•
FPS-Tech	•	•
Leak Reminder	•	•
Altitude Adjustment	•	•
Leak Compensation	•	•
Data Management	•	•
Wireless Transfer	•	•
Cloud Platform	•	•
G.W./N.W.	2.9/1.5kg	
Packing	1SET/CTN 41.5×30×16cm	

BREATHCARE II

AUTO CPAP / BI-LEVEL



NEW AIRWAY / ENJOY BREATHING

◦ New Airway Design

A new airway design that complies with fluid mechanics to improve air duct performance. It increases the flow rate of the whole device by 20% and reduces the impedance by 10%, ensuring the user more air inhaled and an easier process to exhale, thus improving the comfort level.

◦ Magnetic Field Orientation Control System

The upgraded motor drive control system allows no-break speed shifting, quickly responds to the inhale pressure, and enables real-time air leakage compensation.

◦ Acoustic Noise Reduction Design

The sound and audio adjustment technology is used to adjust the high frequency noise to the sound frequency that is inaudible to the human ear. At the same time, the motor is wrapped around the shock absorber structure, enabling multiple acoustic materials to absorb the vibration to further reduce the noise.

YUWELL BREATHCARE PAP / PRODUCT PARAMETERS

Specification	BreathCare II				
Picture					
	White		Black		
Color	White	Black	Black	White	White
Model	YH-450	YH-480	YH-820	YH-825	YH-830
Certification	CE	CE	CE	CE	CE
Mode	CPAP、Auto CPAP	CPAP、Auto CPAP	CPAP、S、T、ST、VGPS	CPAP、S、T、ST、VGPS	CPAP、S、T、ST、VGPS
Pressure	4-20cmH ₂ O	4-20cmH ₂ O	4-20cmH ₂ O	4-25cmH ₂ O	4-30cmH ₂ O
Noise	28dBA±2dBA	28dBA±2dBA	28dBA±2dBA	28dBA±2dBA	28dBA±2dBA
Efficient Humidifier	•	•	•	•	•
Smart Humidification	•	•	•	•	•
Smart Start/Stop	•	•			
FPS-Tech	•	•			
Central Apnea Detection	•	•	•	•	•
Leak Reminder	•	•	•	•	•
High Pressure Reminder			•	•	•
Asphyxiation Reminder			•	•	•
Trigger Adjustment			•	•	•
Cycle Adjustment			•	•	•
Slope Adjustment			•	•	•
Tidal Volume Adjustment			•	•	•
Altitude Adjustment	•	•	•	•	•
Leak Compensation	•	•	•	•	•
Data Management	•	•	•	•	•
Wireless Transfer	•	•	•	•	•
Cloud Platform	•	•	•	•	•
G.W./N.W.	3.6/2.4kg				
Packing	1SET/CTN 40x27x17.5cm				

BREATHCARE | CPAP / AUTO CPAP / BI-LEVEL



STABLE AND RELIABLE QUALITY

- High Performance Motor
 - Silent DC brushless motor
 - 0.2 seconds quick response
 - 20000 hours of service life
- No-Spoiler Air Duct Design
 - Accurately monitors the air flow signals and reduces noise at the same time
- Patented Micro Airflow Laminar Flow Tube
 - Equipped with the world's advanced pressure and flow sensors
 - Accurately collects pressure & flow signal changes
 - Real-time control
- Field-Oriented Control (FOC) Control Algorithm
 - The back propagation pressure control algorithm enables real-time air leakage compensation

YUWELL BREATHCARE PAP / PRODUCT PARAMETERS

Specification	BreathCare I				
Picture	 White&Gray		 Black&Silver		
Color	White&Gray	White&Gray Black&Silver	White&Gray	White&Gray	White&Gray
Model	YH-360	YH-560/580	YH-720	YH-725	YH-730
Certification	CE	CE	CE	CE	CE
Mode	CPAP	CPAP, Auto CPAP	CPAP, S, T, ST, VGPS	CPAP, S, T, ST, VGPS	CPAP, S, T, ST, VGPS
Pressure	4-20cmH ₂ O	4-20cmH ₂ O	4-20cmH ₂ O	4-25cmH ₂ O	4-30cmH ₂ O
Noise	28dBA±2dBA	28dBA±2dBA	28dBA±2dBA	28dBA±2dBA	28dBA±2dBA
Efficient Humidifier	•	•	•	•	•
Smart Start/Stop	•	•			
FPS-Tech	•	•			
Leak Reminder	•	•	•	•	•
High Pressure Reminder			•	•	•
Asphyxiation Reminder			•	•	•
Trigger Adjustment			•	•	•
Cycle Adjustment			•	•	•
Slope Adjustment			•	•	•
Tidal Volume Adjustment			•	•	•
Altitude Adjustment	•	•	•	•	•
Leak Compensation	•	•	•	•	•
Data Management	•	•	•	•	•
G.W./N.W.	3.6/2.4kg				
Packing	1SET/CTN 40x27x17.5cm				

HFNC HF-75A

(HIGH FLOW NASAL CANNULA OXYGEN THERAPY)

- 

Humidifier dry alarm
- 

Oxygen source pressure High / Low alarm
- 

Tube falling alarm
- 

FiO₂ High / Low alarm
- 

Room temp. High / Low alarm
- 

Free of disinfection

- Mode: High / Low Flow
- Flow range: 2-75LPM
- Temperature: 34°C / 31°C-37°C
- FiO₂: 21%-100%
- Screen size: 7-inch touch screen
- Real-time monitoring: Flow, Temperature, FiO₂
- Trend chart: 7 days



Specification	HFNC	
Model	HF-70A	HF-75A
Mode	High flow / Low flow	High flow / Low flow
Flow range	2-25LPM/10-70LPM	2-25LPM/10-75LPM
Temperature	34°C/31°C、34°C、37°C	34°C/31°C、32°C、33°C、34°C、35°C、36°C、37°C
FiO ₂	21-100%	21-100%
SpO ₂ sensor	Optional	Optional
Screen size	7-inch touch screen	7-inch touch screen
Operating mode	Touch	Touch
Disinfection	Free of disinfection	Free of disinfection
Realtime monitoring	Flow, Temperature, FiO ₂ , SpO ₂ (Optional) , Pulse rate(Optional)	Flow, Temperature, FiO ₂ , SpO ₂ (Optional) , Pulse rate(Optional)
Alarm	Humidifier dry alarm Tube falling off alarm Room temperature High / Low alarm Oxygen source pressure High / Low alarm FiO ₂ High / Low alarm Blockage alarm	Humidifier dry alarm Tube falling off alarm Room temperature High / Low alarm Oxygen source pressure High / Low alarm FiO ₂ High / Low alarm Blockage alarm

HFNC HF-60M

(HIGH FLOW NASAL CANNULA OXYGEN THERAPY)

- 

Humidifier dry alarm
- 

Use with oxygen concentrator/cylinder
- 

Tube falling alarm
- 

FiO₂ High / Low alarm
- 

High / Low ambient temperature alarm
- 

Free of disinfection

- Mode: High / Low Flow
- Flow range: 2-60LPM
- Temperature: 34°C / 31°C, 34°C, 37°C
- FiO₂: 21%-100%
- Screen size: 7-inch touch screen
- Heated Humidifier: Reusable and disinfectable
- Real-time monitoring: Flow, Temperature, FiO₂
- Trend chart: 7 days



Specification	HFNC	
Model	HF-60M	HF-70M
Mode	High flow / Low flow	High flow / Low flow
Flow range	2-25LPM/10-60LPM	2-25LPM/10-70LPM
Temperature	34°C/31°C、34°C、37°C	34°C/31°C、34°C、37°C
FiO ₂	21-100%	21-100%
SpO ₂ sensor	Optional	Optional
Screen size	7-inch touch screen	7-inch touch screen
Operating mode	Touch	Touch
Disinfection	Free of disinfection	Free of disinfection
Realtime monitoring	Flow, Temperature, FiO ₂ , SpO ₂ (Optional) , Pulse rate(Optional)	Flow, Temperature, FiO ₂ , SpO ₂ (Optional) , Pulse rate(Optional)
Alarm	Humidifier dry alarm Tube falling off alarm Room temperature High / Low alarm FiO ₂ High / Low alarm Blockage alarm	Humidifier dry alarm Tube falling off alarm Room temperature High / Low alarm FiO ₂ High / Low alarm Blockage alarm

YUWELL BREATHWEAR MASK SERIES



Fit



Comfort



Convenience



Light & Soft




Rebound Silicone Cushion

The mask is made of rebound silicone cushion, which is airtight, light and soft. It is designed to reduce facial pressure to make it fit and comfortable.



Adaptive Silicone Cushion

The silicone cushion has an adaptive design at the nose bridge, which can be automatically adjusted according to different face shapes, so as to better relieve the pressure on the high nose and improve user experience.



Quiet Vent-Hole

Honeycomb design to reduce airflow erosion intensity. Gentle and quiet.



Quick-Release Connector

Using a quick release connector to easily separate or connect the mask and the device.



Soft and Light Headgear

Selected soft skin-friendly fabrics, light and stable, help patients adapt to the treatment comfortably.



TECHNICAL SPECIFICATIONS				
Picture				
Model	YF-01	YN-02	YF-03	YN-01
Certification	CE&FDA	CE&FDA	CE&FDA	CE&FDA
Type	Full Face Mask	Nasal Mask	Full Face Mask	Nasal Mask
Therapy Pressure Range	4-30 cmH ₂ O	4-30 cmH ₂ O	4-30 cmH ₂ O	4-30 cmH ₂ O
Size	S/M/L/XL	S/M/L	S/M/L	S/M/L
Tube Connection	φ22mm	φ22mm	φ22mm	φ22mm
Flow Resistance(At 50 L/min)	≤0.5cmH ₂ O	≤0.5cmH ₂ O	≤0.5cmH ₂ O	≤0.5cmH ₂ O
Flow Resistance(At 100 L/min)	≤1.5cmH ₂ O	≤1.5cmH ₂ O	≤1.5cmH ₂ O	≤1.5cmH ₂ O
Expiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min)	≤2.0cmH ₂ O/(L/S)	—	≤2.0cmH ₂ O/(L/S)	—
Inspiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min)	≤2.0cmH ₂ O/(L/S)	—	≤2.0cmH ₂ O/(L/S)	—
Anti Asphyxia Valve Open Pressure	≤2.0cmH ₂ O/(L/S)	—	≤2.0cmH ₂ O/(L/S)	—
Anti Asphyxia Valve Closed Pressure	≤2.0cmH ₂ O/(L/S)	—	≤2.0cmH ₂ O/(L/S)	—
G.W./N.W.	147/106g	90/60g	139/103g	101/63g
Packing	10PCS/CTN 35×28×28cm	10PCS/CTN 35×24×24cm	10PCS/CTN 35×28×28cm	10PCS/CTN 35×24×24cm

TECHNICAL SPECIFICATIONS			
Picture			
Model	YF-02	YN-03	YP-01
Certification	CE&FDA	CE&FDA	CE&FDA
Type	Full Face Mask	Nasal Mask	Pillow Mask
Therapy Pressure Range	4-30 cmH ₂ O	4-30 cmH ₂ O	4-20 cmH ₂ O
Size	S/M/L/XL	S/M/L	S/M/L
Tube Connection	φ22mm	φ22mm	φ22mm
Flow Resistance(At 50 L/min)	≤0.5cmH ₂ O	≤0.5cmH ₂ O	≤0.5cmH ₂ O
Flow Resistance(At 100 L/min)	≤1.5cmH ₂ O	≤1.5cmH ₂ O	≤1.5cmH ₂ O
Expiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min)	≤2.0cmH ₂ O/(L/S)	—	—
Inspiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min)	≤2.0cmH ₂ O/(L/S)	—	—
Anti Asphyxia Valve Open Pressure	≤2.0cmH ₂ O/(L/S)	—	—
Anti Asphyxia Valve Closed Pressure	≤2.0cmH ₂ O/(L/S)	—	—
G.W./N.W.	148/107g	129/99g	268/81g
Packing	10PCS/CTN 35×28×28cm	10PCS/CTN 35×28×28cm	10PCS/CTN 34×22×24.5cm

yuwell

OXYGEN
CONCENTRATORS

RESPIRATORY



PORTABLE OXYGEN CONCENTRATOR

SPIRIT -3



RESPIRATORY


- 0.02s
Respond in 0.02s
- Portable
- 4-grade regulating
- Detachable molecular sieve
- German technical compressor
- Generate oxygen by Impulses
- Modular lithium battery
- Automatic warning system

- Maximum Flow: 750ml/min (Equivalent Continuous Oxygen Supply: 3.75L/min)
- Oxygen Concentration (Measured after warm up 30minutes): ≥90%
- Maximum Outlet Pressure: 150kPa
- Sound Level: ≤55dB
- Minimum Operating Time: 30minutes
- Operational temperature: 5°C ~ 40°C
- Altitude: 0 to 5000m

QTN: 1PCS/CTN
G.W./N.W: 4.5kgs/2.2kg
Size: 32×32×21cm



Applicable scene



Indoor



Outdoor



Inside the car



Altitude plateau

9F-5B/9F-5BW

- Ultra-silence
- Medical grade
- High oxygen concentration
- Castors
- Fashionable appearance
- Failure detection
- Nebulizer (9F-5BW only)

- Flow range: 0.5~5L/min
- Oxygen concentration: 95.5%-87%
- Maximum outlet pressure: 70kPa
- Sound level: 42dB(A)
- Nebulizing rate: ≥0.1mL/min(9F-5BW only)
- Input power: 400VA
- Minimum operating time: 30minutes

QTN: 1PCS/CTN
G.W./N.W: 20.5/17.5kg
Size: 48×40×63cm

Ultra silence



8F-5A



Global sales leader



Medical grade



High oxygen concentration



5L oxygen flow



Failure detection



Ultra-silence



Castors

- Flow Range: 0.5~5L /min
- Oxygen Concentration: 95.5%-87%
- Maximum Limited Pressure: 70kPa
- Sound Level: 48 dB(A)
- Input Power: 350VA
- Minimum Operating Time: 30minutes

QTN: 1PCS/CTN
G.W./N.W: 17.5/15.5kg
Size: 48.5×34×60cm



Oxygen Concentrator



8F-10(CE)/8F-10W Single flow

- Medical grade
- High oxygen concentration
- 10L oxygen flow
- Cumulative usage time
- Failure detection
- Castors
- HEPA filter
- Nebulizer (8F-10W only)

- Nebulizing function: (8F-10W)
- Maximum Nebulizing Rate: ≥0.15mL/min(8F-10W)
- Automatic reminder of filter cleaning: 100 hours
- Flow range: 1-10L/min
- Oxygen concentration: 95.5%-87%
- Maximum limited pressure: 100kPa
- Sound level: 55dB(A)
- Input power: 700VA
- Minimum operating time: 30minutes

QTN: 1PCS/CTN
N.W: 23.5kg
Size: 41×37.5×66.8cm



8F-10/8F-10W Dual flow

- Medical grade
- High oxygen concentration
- 10L oxygen flow
- Cumulative usage time
- Failure detection
- Castors
- HEPA filter
- Nebulizer (8F-10W only)

- Nebulizing function: (8F-10W)
- Maximum Nebulizing Rate: ≥0.15mL/min(8F-10W)
- Dual flow outlets for two persons
- Automatic reminder of filter cleaning: 100 hours
- Flow range: 1-10L/min (The total flow rate is 10L at dual flow)
- Oxygen concentration: 95.5%-87%
- Maximum limited pressure: 100kPa
- Sound level: 55dB(A)
- Input power: 700VA
- Minimum operating time: 30minutes

QTN: 1PCS/CTN
N.W: 23.5kg
Size: 41×37.5×66.8cm



7F-5/7F-10(CE)

- Failure detection
- Medical grade
- High oxygen concentration
- Castors
- Nebulizer (7F-5W only)

- Flow Range: 0.5~ 5L /min(7F-5) 1~ 10L /min(7F-10)
- Oxygen Concentration: 95.5%-87%
- Maximum Limited Pressure: 70kPa
- Sound Level: 55 dB (A)(7F-5)
≤60dB(A)(7F-10)
- Input Power: 500 VA(7F-5) 850VA(7F-10)
- Minimum Operating Time: 30minutes

QTN: 1PCS/CTN
G.W./N.W: 29/27kg(7F-5)
37.5/33kg(7F-10)
Size: 49.5×43×78cm(7F-5)
57×49.5×80.5cm(7F-10)

Oxygen Concentrator



7F-5B

- Failure detection
- Medical grade
- High oxygen concentration
- Castors
- Nebulizer (7F-5BW only)

- Flow Range: 0.5~ 5L /min
- Oxygen Concentration: 95.5%-87%
- Maximum Limited Pressure: 70kPa
- Sound Level: ≤55 dB (A)
- Input Power: 400 VA
- Minimum Operating Time: 30minutes

QTN: 1PCS/CTN
G.W./N.W: 26/24kg
Size: 56×39.5×71.5cm

Oxygen Concentrator



7F-3NW

- 

Failure detection
- 

Medical grade
- 

High oxygen concentration
- 

Castors
- 

Ultra-silence
- 

Infrared remote control system
- 

Fashionable appearance

- Flow Range: 0.5~3L/min.
- Oxygen Concentration: ≥90%.
- Maximum Limited Pressure: 50kPa
- Sound Level: ≤50dB(A)
- Input Power: 230VA

QTN: 1PCS/CTN
G.W./N.W: 11/9.8kg
Size: 42×26×55cm



7F-3/7F-3W

- 

Failure detection
- 

Medical grade
- 

High oxygen concentration
- 

Castors
- 

Nebulizer
(7F-3W only)

- Flow Range: 0.5~ 3L /min
- Oxygen Concentration: 93%±3%
- Maximum Limited Pressure: 50kPa
- Sound Level: ≤55 dB (A)
- Input Power: 320 VA
- Minimum Operating Time: 30minutes

QTN: 1PCS/CTN
G.W./N.W: 26/24kg
Size: 56×39.5×71.5cm



YU300

- 

Fashionable appearance
- 

Ultra-silence

- Flow Range: 1~ 5L /min
- Oxygen Concentration: 90%~30%
- Sound Level: 43 dB (A)(Front side)
49 dB (A)(Whole unit)
- Input Power: 135 VA
- Minimum Operating Time: 30minutes

QTN: 1PCS/CTN
G.W./N.W: 9.5/7.5kg
Size: 36.5×24.5×38cm

Oxygen Concentrator



8F-1A/8F-2A

- 

Failure detection
- 

Medical grade
- 

High oxygen concentration
- 

Fashionable appearance

- Maximum Recommended Flow: 8F-1A: 1L/min
8F-2A: 2L/min
- Oxygen Concentration: ≥90%
- Maximum Limited Pressure: 100kPa
- Sound Level: ≤50dB(A)
- Input Power: 150VA
- Minimum Operating Time: 30minutes

QTN: 1PCS/CTN
G.W./N.W: 11/9kg
Size: 41×29.5×47cm

yuwell

NEBULIZERS

RESPIRATORY



AIR-COMPRESSING NEBULIZER


403S



- Medical grade
- Nebulizer rate adjustable
- Portable
- Oil-free pump
- Small nebulizer particles
- Easy & visualized operation

- Power supply: ☐AC120V±10% ☐AC220V±10% ☐AC230V±10%
☐50Hz ☐60Hz
 - Max pressure: ≥0.15MPa
 - Max air flow: ≥10L/min
 - Max nebulizing rate: ≥0.2mL/min
 - MMAD: approx.2.5µm, proportion of <5µm: >65%
 - Residual Volume: ≤1mL
 - Operation noise: approx.60dB(A)
 - Operation time: continuous operation
- QTN: 10PCS/CTN
G.W./N.W.: 14.1/11.2kg
Size: 60×47×29.5 cm

Introduction of Accessories



NEW



405K Air-Compressing Nebulizer

- Dual DC pump
- Portable
- Timing
- Easy & visualized operation
- Three levels adjustable flow
- Soundless

- Power supply: AC100V-240V , 50/60Hz
- Max nebulizing rate: ≥0.2mL/min
- Residual Volume: ≤1mL
- Volume of medicine cup: 8ml
- MMD: 3.9um±25%
- Max air flow: ≥3.5L/min
- Operation noise: approx.50dB(A)
- Ingress Protection: IP21
- Operation time: 30 minutes on 40 minutes off
- N.W: approximately 0.67kg
- Dimension:143×152×83mm

QTN: 6PCS/CTN
G.W./N.W.: 7.5/4kg
Size: 60×40.5×22cm

Air-Compressing Nebulizer

NEW



405M Air-Compressing Nebulizer

- Dual DC pump
- Portable
- Timing
- Easy & visualized operation
- Three levels adjustable flow
- Soundless

- Power supply: AC100V-240V, 50/60Hz
- Max nebulizing rate: 0.2ml/min
- Residual Volume: ≤1mL
- Volume of medicine cup: 8ml
- MMD: 3.9um±25%
- Max air flow: ≥3.5L/min
- Operation noise: approx.55dB(A)
- Ingress Protection: IP21
- Operation time: 20 minutes on 40 minutes off
- N.W.: approximately 0.62kg
- Dimension: 150×156×74mm

QTN: 6PCS/CTN
G.W./N.W.: 8.5/3.7kg
Size: 63×42.5×21cm

Air-Compressing Nebulizer

NEW



407C
Air-Compressing Nebulizer

- 
Medical grade
- 
Hospital prefer
- 
Nebulizer rate adjustable
- 
Timing
- 
Nebulizer particle diameter adjustable
- 
Oil-free pump
- 
Non-stop operation
- 
Filter replacement prompt

- Power supply: ☐AC120V±10% ☐AC220V±10% ☐AC230V±10%
☐50Hz ☐60Hz
- Max nebulizing rate: ≥0.2ml/min
- Residual Volume: ≤1mL
- Max air flow: ≥7L/min
- I gear(MMD): 2.8μm±25%
- II gear(MMD): 3.9μm±25%
- III gear(MMD): 5.7μm±25%
- Residual Volume: ≤1.0ml
- Operation noise: ≤60dB (A)
- Ingress Protection: IP21
- Operation time: continuous operation
- N.W.: 2.5kg
- Dimension: 237×186×121mm

QTN: 4PCS/CTN
G.W./N.W.: 14.2/9.6kg
Size: 65.5×49×26cm

CE



403T (CE/Basic sets)
Air-Compressing Nebulizer

- 
Medical grade
- 
Nebulizer rate adjustable
- 
Oil-free pump
- 
Small nebulizer particles
- 
Portable

- Power supply: ☐AC120V±10% ☐AC220V±10% ☐AC230V±10%
☐50Hz ☐60Hz
- Max pressure: ≥0.15mpa
- Max air flow: ≥10L/min
- Max nebulizing rate: ≥0.2mL/min
- MMAD: approx.2.5μm, proportion of <5μm: >65%
- Residual Volume: ≤1mL
- Operation noise: ≤60dB(A) (approx.55dB(A))
- Operation time: continuous operation (CE)
- Operation time: 30 minutes on,30 minutes off (Basic sets)

QTN: 10PCS/CTN
G.W./N.W.: 19/14.7kg(CE)
16.5/11.5kg(Basics sets)
Size: 85×30×40cm

CE



403A(CE)/403AI(Basic sets)
Air-Compressing Nebulizer

- 
Medical grade
- 
Nebulizer rate adjustable
- 
Oil-free pump
- 
Small nebulizer particles
- 
Storage box

- Power supply: ☐AC120V±10% ☐AC220V±10% ☐AC230V±10%
☐50Hz ☐60Hz
- Max pressure: ≥ 0.15MPa
- Max air flow: ≥10L/min
- Max nebulizing rate: ≥0.2mL/min
- MMAD: approx.2.5μm, proportion of <5μm: >65%
- Residual Volume: ≤1mL
- Operation noise: ≤65dB(A)
- Operation time: continuous operation

QTN: 6PCS/CTN
G.W.: 17kg
Size: 70×44×32cm

Air-Compressing Nebulizer



403D
Air-Compressing Nebulizer

- 
Medical grade
- 
Nebulizer rate adjustable
- 
Oil-free pump
- 
Small nebulizer particles
- 
Portable


- Power supply: ☐AC120V±10% ☐AC220V±10% ☐AC230V±10%
☐50Hz ☐60Hz
- Max pressure: ≥0.15MPa
- Max air flow: ≥10L/min
- Max nebulizing rate: ≥0.2mL/min
- MMD: 3.9um±25%
- Residual Volume: ≤1mL
- Operation noise: ≤65dB(A)
- Operation time: 30 minutes on,30 minutes off

QTN: 6PCS/CTN
G.W.: 12.4kg
Size: 70×28.5×32cm

Air-Compressing Nebulizer



405A Air-Compressing Nebulizer




- Input: AC100 -240V, 50/60Hz
- Output: DC12V, 1000mA
- Max nebulizing rate: $\geq 0.2\text{ml/min}$
- MMAD: $3.7\mu\text{m}$, proportion of $<5\mu\text{m}$: $>60\%$
- Max pressure: $\geq 0.1\text{MPa}$
- Max air flow: $\geq 3.5\text{L/min}$
- Residual Volume: $\leq 1\text{mL}$
- Operation noise: approx. 45dB(A)
- Operation time: 20minutes on, 40 minutes off

QTN: 10PCS/CTN
G.W./N.W.: 7.2/6.6kg
Size: 56.5×33×14.5cm



405B Air-Compressing Nebulizer




- Input: AC100 -240V, 50/60Hz
- Output: DC12V, 1000mA
- Max nebulizing rate: $\geq 0.2\text{ml/min}$
- MMAD: $3.7\mu\text{m}$, proportion of $<5\mu\text{m}$: $>60\%$
- Max pressure: $\geq 0.1\text{MPa}$
- Max air flow: $\geq 3.5\text{L/min}$
- Residual Volume: $\leq 1\text{mL}$
- Operation noise: approx 45dB(A)
- Operation time: 20minutes on, 40 minutes off

QTN: 10PCS/CTN
G.W./N.W.: 5.3/2.2kg
Size: 47×32×17.5cm

Air-Compressing Nebulizer



407A Air-Compressing Nebulizer




- Power supply: AC220±22V, 50±1Hz
- Input power: 200VA
- Max pressure: $\geq 0.15\text{MPa}$
- Max air flow: $\geq 10\text{L/min}$
- Max nebulizing rate: $\geq 0.2\text{mL/min}$ (Actual measurement is about 0.4ml/min)
- MMAD: $3.0\mu\text{m}$, proportion of $<5\mu\text{m}$: approx.60%
- Residual Volume: $\leq 1\text{mL}$
- Operation noise: $\leq 60\text{dB(A)}$
- Operation time: continuous operation

QTN: 2PCS/CTN
G.W.: 6.1kg
Size: 54.5×27.5×23.5cm

Ultrasonic Nebulizer



402AI Ultrasonic Nebulizer



- Power supply: AC110~220V, 50/60Hz
- Ultrasonic frequency: $1.7\text{MHz} \pm 10\%$
- Input power: 50VA
- MMD: $3.9\mu\text{m}$
- Maximum nebulizing rate: $\geq 3\text{mL/min}$
- Timing range: 0~60min
- Continuous loading time: $\geq 4\text{h}$
- Fill volume of medicine cup: 150mL
- Fill volume of water container: 350mL

QTN: 2PCS/CTN
G.W./N.W.: 4/3.5kg
Size: 41.5×30×26.5cm

MESH
NEBULIZER
M102(CE) / M118



Medical grade

Portable

Soundless

Small nebulizer particles

- Power supply:
DC 3.0V (two "AA" 1.5V alkaline batteries)(M102)
DC 3.7V(-0.2V, +0.5V) (Build-in lithium battery)(M118)
Input: AC100V~240V 50Hz/60Hz 0.2A;
Output: DC5V 1000 mA
 - Nebulization rate: ≥0.2 mL/min
 - Medication residue: ≤0.5mL
 - Particle Size: MMD 4.2μm±25%
 - Percent of MMD < 5μm: > 50%
 - Noise: ≤45dB(A)
- QTN: 10PCS/CTN
G.W.: 4.7kg
Size: 40.5×33.5×17cm

Introduction of Accessories



M103(CE) / M119
Mesh Nebulizer



Medical grade

Portable

Soundless

Small nebulizer particles

- Input power:
DC 3.0V (two "AA" 1.5V alkaline batteries)(M103)
DC 3.7V(-0.2V, +0.5V) (Build-in lithium battery)(M119)
- Operation noise: approx.40dB(A)
- MMAD: 3.9μm, proportion of <5μm: >50%
- Max nebulizing rate: ≥0.2ml/min
- Residual Volume: ≤0.25mL

QTN: 10PCS/CTN
G.W./N.W.: 4.7/0.97kg
Size: 40.5×33.5×17cm

Mesh Nebulizer



M105
Mesh Nebulizer

Medical grade

Nebulizer rate adjustable

2-grade regulating

Soundless

Portable

Small nebulizer particles

- Input power: 3V(two "AA" 1.5V alkaline batteries).
Input: ~100-240V 50/60Hz 0.2A;
Output: 5V 1000 mA
- Nebulization rate: ≥0.2 mL/min
- Medication residue: ≤0.3 mL
- Particle Size: MMD Approx. 3.7μm
- Noise: ≤45dB(A)


QTN: 10PCS/CTN
G.W./N.W.: 4.7/0.97kg
Size: 40.5×33.5×17cm

yuwell

BLOOD GLUCOSE METERS

DIABETES
MANAGEMENT





660



Auto-coding



Memory capacity



FAD-GDH

QTN: 20PCS/CTN
G.W./N.W.: 5.2/1.0kg
Size: 59x24x18.5cm
Weight: 50g

- Enzyme: FAD-GDH
- Unit: mg/dL or mmol/L
- Measurement range: 0.6-33.3mmol/L(10-600mg/dL)
- Blood volume: 0.8μL
- Hematocrit: 10%-70%
- Measuring time: 8 seconds
- Memory capacity: 500 groups
- Average results: 7days,14days, 30days, 60days and 90days
- Operation temperature: 5℃-40℃
- Display screen: LCD
- Target bloods: Fingertip fresh capillary whole blood or venous whole blood
- Equipment operation mode: Continuous operation
- Input power: DC 3.0V 100mA
- Battery life: 2 AAA alkaline batteries can measure about 1000 times
- Service life: 5 years(Except for wearing parts)

650



Auto-coding



Memory capacity



FAD-GDH

Optional function



Bluetooth


QTN: 20PCS/CTN
G.W./N.W.: 5.2/1.0kg
Size: 59x24x18.5cm
Weight: 50g

- Enzyme: FAD-GDH
- Unit: mg/dL or mmol/L
- Measurement range: 0.6-33.3mmol/L(10-600mg/dL)
- Blood volume: 0.8μL
- Hematocrit: 10%-70%
- Measuring time: 8 seconds
- Memory capacity: 800 groups
- Average results: 7days,14days, 30days, 60days and 90days
- Operation temperature: 5℃-40℃
- Display screen: LCD
- Target bloods: Fingertip fresh capillary whole blood or venous whole blood
- Equipment operation mode: Continuous operation
- Battery life: 2 AAA alkaline batteries can measure about 1000 times
- Service life: 5 years(Except for wearing parts)

NEW



305



Optional function



Bluetooth


- Enzyme: GOD
- Strip reject button
- Unit: mg/dL or mmol/L
- Blood volume: 1 μ L
- Measuring time: 8 seconds
- Measuring range: 0.6-33.3mmol/L(10-600mg/dL)
- Memory capacity: 500 groups
- Average results: 7days,14days, 30days, 60days and 90days
- Operation temperature: 5°C-40°C
- Hematocrit range: 30%-60%
- 2 AAA batteries
- Bluetooth(Optional)

QTN: 20PCS/CTN
G.W./N.W.: 5.25/1.04kg
Size: 59×24×18.5cm
Weight: 50g

NEW



720



Optional function



Bluetooth

- Enzyme: GOD
- Strip reject button
- Unit: mg/dL or mmol/L
- Blood volume: 1 μ L
- Measuring time: 8 seconds
- Measuring range: 0.6-33.3mmol/L(10-600mg/dL)
- Memory capacity: 800 groups
- Average results: 7days,14days, 30days, 60days and 90days
- Operation temperature: 5°C-40°C
- Hematocrit range 30%-60%
- 2 AAA batteries
- Bluetooth(Optional)

QTN: 20PCS/CTN
G.W./N.W.: 5.25/1.04kg
Size: 59×24×18.5cm
Weight: 50g

CE




710



- Enzyme: GOD
- Unit: mg/dL or mmol/L
- Blood volume: 1 μ L
- Measuring time: 8 seconds
- Measuring range: 0.6-33.3mmol/L(10-600mg/dL)
- Memory capacity: 250 test results
- Average results: 7days,14days, 30days, 60days and 90days
- Hematocrit Range: 30%-60%
- 1 Button cell

QTN: 20PCS/CTN
G.W./N.W.: 3.8/0.98kg
Size: 59×24×18.5cm
Weight: 40g

750



Optional function




Bluetooth

- Enzyme: GOD
- Unit: mg/dL or mmol/L
- Blood volume: 1 μ L
- Measuring time: 8 seconds
- Measuring range: 0.6-33.3mmol/L(10-600mg/dL)
- Memory capacity: 500 test results
- Average results: 7days,14days, 30days, 60days and 90days
- Hematocrit Range: 30%-60%
- 2 AAA batteries


QTN: 20PCS/CTN
G.W./N.W.: 4.4/0.8kg
Weight: 49g



582



Optional function
(Choose one of the three)




- Enzyme: GOD
- Strip reject button
- Unit: mg/dL or mmol/L
- Blood volume: 1µL
- Measuring time: 8 seconds
- Measuring range: 0.6-33.3mmol/L(10-600mg/dL)
- Memory capacity: 250 groups
- Average results: 7days,14days, 30days, 60days and 90days
- Operation temperature: 5°C-40°C
- Hematocrit range 30%-60%
- 2 AAA batteries

QTN: 20PCS/CTN
G.W./N.W.: 5.25/1.04kg
Size: 59×24×18.5cm
Weight: 50g



580



Optional function



Backlit LCD display

- Enzyme: GOD
- Metering unit: mg/dL or mmol/L
- Blood volume: 1ul
- Measuring time: 8seconds
- Measurement range: 0.6-33.3mmol/L(10-600mg/dL)
- Memory capacity: 250 groups
- Average results: 7days,14days, 30days, 60days and 90days
- 2 AAA batteries
- Before/after meal function

QTN: 20PCS/CTN
G.W./N.W.: 4.57/1.5kg
Size: 59x24x18.5cm
Weight: 75g



Model	660	650	305	720	750
Voice broadcast					
Backlit LCD display					
Auto-coding	•	•	•	•	•
Average results	•	•	•	•	•
Bluetooth		Optional	Optional	Optional	Optional
Strip eject button	•	•	•	•	•
Enzyme	GDH	GDH	GOD	GOD	GOD
Blood volume	0.8ul	0.8ul	1ul	1ul	1ul
Measuring time	8s	8s	8s	8s	8s
Memory capacity	500 groups	800 groups	500 groups	800 groups	500 groups
Size of the whole case	QTN: 20PCS/CTN G.W./N.W.: 5.2/1kg Size: 59×24×18.5cm Weight: 50g	QTN: 20PCS/CTN G.W./N.W.: 5.2/1kg Size: 59×24×18.5cm Weight: 50g	QTN: 20PCS/CTN G.W./N.W.: 5.25/1.04kg Size: 59×24×18.5cm Weight: 50g	QTN: 20PCS/CTN G.W./N.W.: 5.25/1.04kg Size: 59×24×18.5cm Weight: 50g	QTN: 20PCS/CTN G.W./N.W.: 4.4/0.8kg Weight: 49g
Certification	CE				



Model	710	582	580		
Voice broadcast		Optional			
Backlit LCD display		Optional	Optional		
Auto-coding	•	•	•		
Average results	•	•	•		
Bluetooth		Optional			
Strip eject button		•	•		
Enzyme	GOD	GOD	GOD		
Blood volume	1µl	1ul	1ul		
Measuring time	8s	8s	8s		
Memory capacity	250 groups	250 groups	250 groups		
Size of the whole case	QTN: 20PCS/CTN G.W./N.W.: 3.8/0.98kg Size: 59×24×18.5cm Weight: 40g	QTN: 20PCS/CTN G.W./N.W.: 5.25/1.04kg Size: 59×24×18.5cm Weight: 50g	QTN: 20PCS/CTN G.W./N.W.: 4.57/1.5kg Size: 59×24×18.5cm Weight: 75g		
Certification	CE	CE			

BLOOD GLUCOSE
AND URIC ACID METER

GU100



Two in one



Memory capacity



Average results



Bluetooth



Backlit LCD display



Voice Broadcast

 **Optional function**
(Choose one of the three)



Specification	Blood Glucose	Uric Acid
Strip model	Y330	YU300
Test range	10 mg/dL to 600 mg/dL(0.6 mmol/L to 33.3 mmol/L)	3.04 mg/dL to 19.96 mg/dL(181 μmol/L to 1188 μmol/L)
Target blood	Fresh capillary whole blood	Fresh capillary whole blood and venous whole blood
Accuracy	<100mg/dL:Within ±15mg/dL ≥100mg/dL:Within ±15%	≤ 5.04 mg/dL:Within±1.01mg/dL ≤ 5.04 mg/dL:Within ±20%
Sample volume	1μL	3 μL
Optional function	Voice broadcast、Backlight、Bluetooth	
Shelf life for strip	Unopened:24 months After opening:3 months	Unopened:18 months After opening:3 months
Measuring time	<10 seconds	<10 seconds
Memory groups	250 groups	250 groups
Atmospheric pressure range	50 kPa to 106 kPa	50 kPa to 106 kPa
Storage Ambient temperature range for meter	-10 °C~ +50 °C	-10 °C~ +50 °C
Storage Ambient temperature range for strip	4°C-30°C	4°C-30°C
Testing Ambient temperature range for strip.	10°C-40°C	10°C-40°C
Hematocrit	30%-60%	30%-60%


Packaging parameters

QTN (PCS/CTN)	Meter Size (mm)	Meter Weight (g)	G.W./N.W. (kg)	Size (cm)
20	109 (L)×76 (W)×23.5 (H)	71 (excluding the batteries)	5.7/1.37	59×24×18.5

TEST STRIPS

Y330

- Enzyme:GOD
- CE Certified test strip quality
- 50 pcs (1 container)
- Required volume of blood: 1 ul
- Code-free
- Easy to use
- Medical device
- Must be used within max.3 months of opening and before the expiry date



Y350

- Enzyme:FAD-GDH
- CE Certified test strip quality
- 50 pcs (1 container)
- Required volume of blood: 0.8 ul
- Code-free
- Easy to use
- Medical device
- Must be used within max.3 months of opening and before the expiry date



Uric acid test strip
YU300

- Required volume of blood: 3 ul
- Easy to use
- Medical device
- Must be used within max.3 months of opening and before the expiry date



yuwell

CONTINUOUS GLUCOSE MONITORING SYSTEMS

DIABETES
MANAGEMENT



MEET THE ANYTIME CT3 SYSTEM



01.AUTO-APPLICATOR & SENSOR
The simple auto-applicator is a one-touch applicator that easily inserts a small built-in sensor just beneath the skin.

03.WIRED CHARGING STAND
For one time charge only takes less than 1.5 hours.

02.TRANSMITTER (2-year shelflife)
The rechargeable transmitter can be used for up to two years. It can support 14 days' wearing.

APPLICABLE SCENE



Why should I choose CGM CT3?

- 

No fingersticks
Without the need for fingersticks and painless to wear
- 

Affordable
Reusable transmitter makes it affordable for all users even without insurance coverage
- 

Accurate
Keep tracking for 14-day accuracy with clinical proven
- 

Feel in control
Reduce hypoglycemia day and night, and lower HbA1c levels

PARAMETER

Sensor shelflife	1 year
Sensor use time	14 days
Transmitter shelflife	2 years
Charge time	<1.5h
Waterproof	IP58
Mard	9.07%
App	Anytime, Follow, Glooko
Bluetooth transmission range	6M
Glucose reading frequency	Every three minutes
Alarms	Predictive and realtime alarm
Power supply	Rechargeable
Storage temperature	2-30℃
Transport temperature	2-45℃
QTN	40PCS/CTN

yuwell

BLOOD PRESSURE MONITORS

DIAGNOSTICS



MEDICAL AUTOMATIC ELECTRONIC
BLOOD PRESSURE MONITOR



YE990



- Large Backlit display
- Prompt for movement
- Body induction
- Printing function

- Animated voice prompt, for automatic measurement
- Voice broadcast function, broadcasting the measured blood pressure in a timely manner
- Printer function, printing the measured blood pressure value
- Supporting multiple transmission modes
- Intelligent pressurization and fixed pressurization
- Dual-mode measurement

Model	YE990
Applicable arm circumference	17~45 cm
Pressure measuring range	0~300mmHg(0.0~40.0kPa)
Pulse measuring range	40~200 time/min
Accuracy of pressure measurement	±2mmHg(±0.267kPa)
Accuracy of pulse measurement	±2%or±2 time/min
Measurement time	About 35s
Data transmission	Serial port, Ethernet, USB, Bluetooth, WIFI
Systolic pressure display	●
Diastolic pressure display	●
Pulse rate display	●
Voice broadcast function	●
Printing function	●
Blood pressure level display function	●
Prompt for movement error	●



DIAGNOSTICS

YE900

- 100K 100,000 times of measurements
- Large backlit LCD
- Portable
- A/B Double measuring methods
- Ultra quiet
- Bluetooth (Optional)
- Anti-collision & waterproof
- Optional S M L cuff size

- Large backlit LCD display
- 120 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for movement error
- Cuff storage
- Portable
- Double measuring methods
- Selectable pressuring
- Anti-collision design
- Dual-mode of automatic and auscultation measurement

QTN: 10PCS/CTN
G.W./N.W.: 20.6/5.9kg
Size: 66.5×40.5×33cm
Weight: 584g



YE630CR

- 100K: 100,000 times of measurements
- Large backlit LCD
- Portable
- Ultra quiet
- Rechargeable
- Anti-collision & waterproof

- Indication for movement error
- Detection for cuff wearing
- Voice function
- Bluetooth transmission
- Lithium battery
- 50 groups of memory

QTN: 10PCS/CTN
G.W./N.W.: 2.54/8.42kg
Size: 66.5×28×22cm
Weight: 257g

Upper Arm Blood Pressure

NEW



YE650AR

- LCD digital display with three colors backlit(Optional)
- Voice broadcast function
- 60 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for movement error
- Average of three measurements
- Bluetooth(Optional)
- Lithium Battery

QTN: 10PCS/CTN
G.W./N.W.: 7.3/3.6kg
Size: 61.5×28.5×22cm
Weight: 270g



YE680B

- Voice broadcast function
- Large backlit LCD display
- 99×2 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for movement error
- Average of three measurements
- Bluetooth(Optional)

QTN: 10PCS/CTN
G.W./N.W.: 9.6/3.3kg
Size: 61.5×28.5×22cm
Weight: 310g

Upper Arm Blood Pressure



YE670A

- Voice broadcast function(Optional)
- (Voice can be turned off)
- Large LCD display
- 74×2 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for movement error
- Average of three measurements

QTN: 10PCS/CTN
G.W./N.W.: 8.3/3.8kg
Size: 50×30×22cm
Weight: 294g



YE670CR

- Voice broadcast function
(Voice can be turned off)
- Large LCD display
- 60×2 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for movement error
- Average of three measurements
- Bluetooth(Optional)
- Lithium Battery

QTN: 10PCS/CTN
G.W./N.W.: 7.6/4.0kg
Size: 50×30×22cm
Weight: 303g



YE660D(CE)/YE660AR

- Large LCD display
- 90 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for movement error
- Average of three measurements
- Lithium Battery(YE660AR only)

QTN: 10PCS/CTN
G.W./N.W.: 8.8/3.5kg
Size: 50×30×22cm
Weight: 260g



YE660E

- Voice broadcast function
- Large LCD digital display
- 99 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for movement error
- Average of three measurements

QTN: 10PCS/CTN
G.W./N.W.: 7.3/3.6kg
Size: 61.5×28.5×22cm
Weight: 270g

Upper Arm Blood Pressure



YE660F

- Voice broadcast function
- Large Backlit LCD display
- 90 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for movement error
- Average of three measurements

QTN: 10PCS/CTN
G.W./N.W.: 8.7/3kg
Size: 61.5×28.5×22cm
Weight: 300g



YE620B

- Voice broadcast function
- Large Backlit LCD display
- 99 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for movement error
- Average of three measurements

QTN: 10PCS/CTN
G.W./N.W.: 8.4/2.9kg
Size: 61.5×28.5×22cm (with 360 cuff)
Weight: 290g



YE610D(CE)/YE610AR

- 80 groups of memory
- Average of three measurements
- Lithium Battery(YE610AR only)

QTN: 10PCS/CTN
G.W./N.W.: 8.8/3.5kg
Size: 50×30×22cm
Weight: 251g



YE610A

- Large LCD digital display(Optional)
- Voice broadcast function
- 50 groups of memory
- Average of three measurements
- Systolic、Diastolic and Pulse rate display

QTN: 10PCS/CTN
G.W./N.W.: 8.8/3.5kg
Size: 50×30×22cm
Weight: 261g

Wrist Blood Pressure



YE8800AR

- Voice broadcast function
- Indication for Position
- Large LCD display
- 74 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for Movement Error
- Average of three measurements
- Lithium Battery
- Bluetooth(Optional)

QTN: 10PCS/CTN
G.W./N.W.: 3.1/1.1kg
Size: 53×24.5×13cm
Weight: 109g



YE8300B

- Large LCD display
- 99 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for Movement Error
- Average of three measurements
- Automatically shutdown in 3 minutes after measurement

QTN: 10PCS/CTN
G.W./N.W.: 2.5/1kg
Size: 36.5×20.5×13cm
Weight: about99g



YE8900A

- Indication for Position
- Large LCD display
- 60 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for Movement Error
- Average of three measurements
- Bluetooth(Optional)

QTN: 10PCS/CTN
G.W./N.W.: 2.4/0.9kg
Size: 38.5×23×11.5cm
Weight: 90g



YE8600A

- Indication for Position
- Large LCD display
- 60 groups of memory
- Indication for irregular heartbeats
- Detection for cuff wearing
- Indication for Movement Error
- Average of three measurements

QTN: 10PCS/CTN
G.W./N.W.: 2.4/0.9kg
Size: 38.5×23×11.5cm
Weight: 90g

Blood Pressure

Upper Arm


Model	YE610A	YE610D	YE660D	YE660E	YE660F
Voice broadcast	Optional		Optional	●	●
Backlit LCD display					●
Bluetooth					
Chargeable					
Indication for irregular heartbeats			●	●	●
Detection for cuff wearing			●	●	●
Indication for Movement Error			●	●	●
Average of three measurements	●	●	●	●	●
Memory groups	50 groups	80 groups	90 groups	99 groups	90 groups
External power supply: Regular	USB Type-C	USB Type-C	USB round hole	USB round hole	Adapter
External power supply: CE		USB Type-C	Adapter(Optional)	Adapter(Optional)	USB round hole
Size of the whole case					
QTN	10PCS/CTN	10PCS/CTN	10PCS/CTN	10PCS/CTN	10PCS/CTN
G.W./N.W.	8.8/3.5kg	8.8/3.5kg	8.8/3.5kg	7.3/3.6kg	8.7/3.0kg
Size	50×30×22cm	50×30×22cm	50×30×22cm	61.5×28.5×22cm	61.5×28.5×22cm
Weight	261g	251g	260g	270g	300g
Certification	CE		CE	CE&FDA	CE&FDA

Wrist


Model	YE8300B	YE8600A	YE8900A	YE8800AR
Voice broadcast				●
Large LCD display	●	●	●	●
Bluetooth			●	●
Lithium battery				●
Indication for position		●	●	●
Irregular heartbeats detection		●	●	●
Detection for cuff wearing		●	●	●
Indication for Movement Error		●	●	●
Average of three measurements	●	●	●	●
Memory groups	99 groups	60 groups	60 groups	74 groups
External power supply	None	None	None	USB(Type C)
Size of the whole case				
QTN	10PCS/CTN	10PCS/CTN	10PCS/CTN	10PCS/CTN
G.W./N.W.	2.5/1.0kg	2.4/0.9kg	2.4/0.9kg	3.1/1.1kg
Size	36.5×20.5×13cm	38.5×23×11.5cm	38.5×23×11.5cm	53×24.5×13cm
Weight	about 99g	90g	90g	109g
Certification	CE			CE&FDA

					
YE620B	YE650AR	YE670A	YE670CR	YE680B	YE900
●	Optional	Optional	●	●	
●	Optional			●	●
	Optional		Optional	Optional	Optional
	●		●		●
●	●	●	●	●	●
●	●	●	●	●	●
●	●	●	●	●	●
●	●	●	●	●	
99 groups	60 groups	74x2 groups	60x2 groups	99x2 groups	120 groups
Adapter	USB Type-C	USB Type-C	USB Type-C	Adapter	Adapter
USB round hole	USB Type-C	USB Type-C	USB Type-C	Adapter(Optional)	
10PCS/CTN	10PCS/CTN	10PCS/CTN	10PCS/CTN	10PCS/CTN	10PCS/CTN
8.4/2.9kg	6.6/4.4kg	8.3/3.8kg	7.6/4.0kg	9.6/3.3kg	20.6/5.9kg
61.5×28.5×22cm	50×30×22cm	50×30×22cm	50×30×22cm	61.5×28.5×22cm	66.5×40.5×33cm
290g	306g	294g	303g	310g	584g
CE&FDA	CE	CE	CE	CE&FDA	


CORRECT MEASUREMENT POSTURE



Remove the thicker clothing from the arm, expose the upper arm or leave only the thinner clothing. Place the cuff on the detected upper arm so that the trachea is placed inside the arm and the cuff nozzle is pointed to the lower arm. Make sure the lower edge of the cuff is 2~3 cm away from the elbow joint.



Place your lower arms flat on the table, palms up, and keep body straight. Note that the rubber airway on the cuff does not twist or break.



- 1. Keep body straight, sit on the chair, the ideal height difference between the table and the chair is 25~30 centimeters.
- 2. Place the center of the cuff at the same height as the heart.

PORTABLE VITAL SIGNS MONITOR

HM100

Measurable parameters:

- ①Blood pressure
- ②SpO2
- ③Body temperature
- ④Blood glucose
- ⑤Uric acid



NEW

Model	HM100
Applicable arm circumference	22~32 cm
Pressure measuring range	0~300mmHg(0.0~40.0kPa)
Pulse measuring range	40~200 time/min
SpO ₂ measuring range	0% ~ 100%
Body temperature measuring range	32.0℃ ~ 43.0℃
Blood glucose measuring range	0.6 mmol/L ~ 33.3 mmol/L (10mg/dL ~ 600mg/dL)
Uric acid measuring range	181 μmol/L ~ 1188 μmol/L
Data transmission	Bluetooth, WIFI 4G
Voice broadcast function	●
2 users	●
Indication for irregular heartbeats	●
Prompt for movement error	●

MEDICAL VITAL SIGNS MONITOR

CM100

Measurable parameters:

- ①Blood pressure
- ②SpO2
- ③Body temperature
- ④Pulse rate



NEW

Model	CM100
Applicable arm circumference	17~47cm
Blood pressure measuring range	0~300mmHg(0.0~40.0kPa)
SpO ₂ measuring range	70%-100%±2%
Pulse rate measuring range	25bpm~250bpm±1bpm
Infrared thermometer temp. range	32.0 °C~43.0 °C
Infrared ear thermometer temp. range	20.0 °C~42.2 °C
Memory groups	5000 groups
Data transmission	Bluetooth, WIFI,4G, USB,RS-232
Measuring time	≤15 seconds
Battery running time	22h
Size	147*190*243mm
Large LCD display	●
PI display	●
Lithium battery	●

yuwell

THERMOMETER

DIAGNOSTICS





YT-3 Infrared Thermometer



- Body measurement range: 32°C - 43°C
- Object measurement range: 0°C - 60°C
- Measure distance: 0-5cm
- Time interval of each measurement: ≤1s
- Automatic shutdown time: <2min
- Measurement duration: ≤1s
- Memory retrieval: 12 memory groups

QTN: 40PCS/CTN
G.W./N.W.: 11/3kg
Size: 62.5x40x16.5cm



YT-1 Infrared Thermometer



- Measurement range: 32-43°C
- Temperature unit: °C/°F
- Measurement duration: 1 second
- Automatic shutdown time: 30 seconds
- Measure distance: 0-5cm

QTN: 40PCS/CTN
G.W./N.W.: 8.8/3kgs
Size: 57.5x36x16.5cm

Infrared Thermometer



YT-2 Infrared Thermometer



- Measurement range: 32-43°C
- Temperature unit: °C/°F
- Measurement duration: 1 second
- Automatic shutdown time: 30 seconds
- Measure distance: 0-5cm
- Memory retrieval: 10 memory groups

QTN: 40PCS/CTN
G.W./N.W.: 8.8/3kgs
Size: 57.5x36x16.5cm

Infrared Thermometer



YT-1C
Infrared Thermometer

- Fast measurement
- Backlit display
- Silent mode selection
- Object temperature measurement

- Measurement range: 32-43°C
- Temperature unit: °C/°F
- Memory retrieval: 9 memory groups
- Measurement duration: 1 second
- Automatic shutdown time: 60 seconds
- Measure distance: 0-5cm
- Measuring Color Indicator

QTN: 40PCS/CTN
G.W./N.W.: 5.9/2.4kg
Size: 57.5×36×16.5cm



YT-2C
Infrared Thermometer

- Fast measurement
- Backlit display
- Object temperature measurement

- Body measurement range: 32°C-43°C
- Object measurement range: 0°C-60°C
- Temperature unit: °C/°F
- Memory retrieval: 30 memory groups
- Measurement duration: 1 second
- Automatic shutdown time: 60 seconds
- Measure distance: 0-5cm

QTN: 40PCS/CTN
G.W./N.W.: 8.8/3kg
Size: 57.5×36×16.5cm



YHW-5(CE) / YHW-6 (Bluetooth)
Infrared Thermometer

- Fast measurement
- Automatic shutdown
- Fever warning
- Backlit display
- Bluetooth (YHW-6 only)

- Body measurement range: 32°C-43°C
- Object measurement range: 0°C-60°C
- Measure distance: 0-3cm
- Temperature unit: °C/°F
- Time interval of each measurement: ≤1s
- Automatic shutdown time: 60s±10s
- Measurement duration: ≤1s
- Memory retrieval: 20 memory groups

QTN: 40PCS/CTN
G.W./N.W.: 8.4/2.8kg
Size: 57.5×35.1×15.4cm



DIAGNOSTICS



YHT100 Infrared Ear Thermometer




- Body measurement range: 34°C -42.2°C
- Temperature unit: °C / °F
- Time interval of each measurement: ≤12s
- Automatic shutdown time: 60s±10s
- Measurement duration: ≤2s
- Number of earmuffs: 20
- Memory retrieval: 10 memory groups

QTN: 40PCS/CTN
G.W./N.W.: 10/3kg
Size: 67×37.5×26.5cm



YHT103 Infrared Ear Thermometer




- Body measurement range: 34°C-42.2°C
- Temperature unit: °C/°F
- Time interval of each measurement: ≤12s
- Automatic shutdown time: 60s±10s
- Measurement duration: ≤2s
- Number of earmuffs: 20
- Memory retrieval: 6 memory groups

QTN: 40PCS/CTN
G.W./N.W.: 10/3kg
Size: 67×37.5×26.5cm

Infrared Ear Thermometer



YHT101 Infrared Ear Thermometer



- Body measurement range: 34°C-42.2°C
- Temperature unit: °C/°F
- Time interval of each measurement: ≤12s
- Automatic shutdown time: 60s±10s
- Measurement duration: ≤2s
- Number of earmuffs: 20
- Memory retrieval: 10 memory groups

QTN: 40PCS/CTN
G.W./N.W.: 7.1/2.5kg
Size: 66×36.5×17.5cm

Infrared Ear Thermometer



YHT200
Infrared Ear Thermometer

- 

Fast measurement
- 

Automatic shutdown
- 

Fever warning
- 

Backlit display
- 

Memory retrieval

- Body measurement range: 34°C-42.2°C
- Temperature unit: °C/°F
- Time interval of each measurement: ≤12s
- Automatic shutdown time: 60s±10s
- Measurement duration: ≤2s
- Number of earmuffs: 20
- Memory retrieval: 6 memory groups

QTN: 40PCS/CTN
G.W./N.W.: 9.2/2.9kg
Size: 66×36.5×17.5cm



DIAGNOSTICS

YT310 / YT311
(Oral use & Armpit use)



YT310

YT311



*Storage cover

- YT310 (normal) YT311 (soft)
- Power supply: internal power supply equipment, alkaline button battery (LR41, -10%~+5%)
- Waterproof rating: IP22
- Battery life: about 1000 times (measured 5-10 minutes/time)
- Measuring range: 32.0°C ~ 42.9°C
- Unit of measurement: °C/°F
- Memory function: memorizes the last measured value
- High temperature reminder: if the measured temperature is higher than 37.8°C, it will send out continuously and quickly "Di-di-di" alarm prompt

QTN: 50PCS/CTN
G.W./N.W.: 4.0/1.1kg
Size: 45×34×16cm

YT312 / YT313
(Oral use & Armpit use)



YT312

YT313



*Storage cover

- YT312 (normal) YT313 (soft)
- Power supply: internal power supply equipment, alkaline button battery (LR41, -10%~+5%)
- Waterproof rating: IP22
- Battery life: about 1000 times (measured 5-10 minutes/time)
- Measuring range: 32.0°C ~ 42.9°C
- Unit of measurement: °C/°F
- Memory function: memorizes the last measured value
- High temperature reminder: if the measured temperature is higher than 37.8°C, it will send out continuously and quickly "Di-di-di" alarm prompt

QTN: 50PCS/CTN
G.W./N.W.: 3.9/1.0kg
Size: 45×34×16cm

yuwell

FINGER PULSE OXIMETER

DIAGNOSTICS





YX102



YX110



Automatic
shutdown



Pulse column
display



Bluetooth



LED backlit
display

- Measurement range of pulse rate: 25bpm~250bpm
- Accuracy of SpO2: 70%~100%: $\pm 2\%$
- Accuracy of pulse rate:
 $\pm 1\%$ / $\pm 1\text{bpm}$ (The result depends on the larger value)
- Power consumption: $\leq 40\text{mA}$
- Operation environment: Temperature($^{\circ}\text{C}$): 5~40
Relative humidity: $\leq 80\%$
- Boot method: Open angle to boot
- Power supply: 2x1.5V AAA alkaline batteries
- Display mode: LED
- Dimension: 60x38x35mm
- Weight: 38g (without batteries)
- Automatic shutdown time: 8s

QTN: 40PCS/CTN
G.W./N.W.: 4.1/1.5kg
Size: 45x26x15cm

YX306



YX310



Automatic
shutdown



OLED
display



Pulse wave
display



Bluetooth



Gravity
sense function



Pulse column
display



PI
display

- Measurement range of pulse rate: 25bpm~250bpm
- Accuracy of SpO2: 70%~100%: $\pm 2\%$
- Accuracy of pulse rate:
 $\pm 1\%$ / $\pm 1\text{bpm}$ (The result depends on the larger value)
- Power consumption: $\leq 40\text{mA}$
- Operation environment: Temperature($^{\circ}\text{C}$): 5~40
Relative humidity: $\leq 80\%$
- Boot method: Push button to boot
- Power supply: 2x1.5V AAA alkaline batteries
- Dimension: 60x38x35mm
- Weight: 38g (without batteries)
- Automatic shutdown time: 8s

QTN: 40PCS/CTN
G.W./N.W.: 4.1/1.5kg
Size: 45x26x15cm



YX301




- Measurement range of pulse rate: 25bpm-250bpm
- Accuracy of SpO₂: 70%-100%: ±2%
- Accuracy of pulse rate:
±1% / ±1bpm (The result depends on the larger value)
- Power consumption: ≤30mA
- Operation environment: Temperature(°C): 5~40
Relative humidity: ≤80%
- Boot method: Push button to boot
- Power supply: 2×1.5V AAA alkaline batteries
- Display mode: OLED
- Dimension: 66×35×30mm
- Weight: 27g (without batteries)
- Automatic shutdown time: 8s

QTN: 40PCS/CTN
G.W./N.W.: 3.8/1.2kg
Size: 45×26×15cm



YX105




- Measurement range of pulse rate: 25bpm-250bpm
- Accuracy of SpO₂: 70%-100%: ±2%
- Accuracy of pulse rate:
±1% / ±1bpm (The result depends on the larger value)
- Power consumption: ≤40mA
- Operation environment: Temperature(°C): 5~40
Relative humidity: ≤80%
- Boot method: Push button to boot
- Power supply: 2×1.5VAAA alkaline batteries
- Display mode: LED
- Dimension: 58×33×33.5mm
- Weight: 35g (without batteries)
- Automatic shutdown time: 8s

QTN: 40PCS/CTN
G.W./N.W.: 4.0/1.3kg
Size: 45×26×15cm



YX106



- Measurement range of pulse rate: 25bpm-250bpm
- Accuracy of SpO₂: 70%-100%: ±2%
- Accuracy of pulse rate:
±1% / ±1bpm (The result depends on the larger value)
- Power consumption: ≤40mA
- Operation environment: Temperature(°C): 5~40
Relative humidity: ≤80%
- Boot method: Push button to boot
- Power supply: 2×1.5VAAA alkaline batteries
- Display mode: LED
- Dimension: 58×33×33.5mm
- Weight: 35g (without batteries)
- Automatic shutdown time: 8s

QTN: 40PCS/CTN
G.W./N.W.: 4.3/1.6kg
Size: 45×26×15cm



YX103



- Measurement range of pulse rate: 25bpm-250bpm
- Accuracy of SpO₂: 70%-100%: ±2%
- Accuracy of pulse rate:
±1% / ±1bpm (The result depends on the larger value)
- Power consumption: ≤40mA
- Operation environment: Temperature(°C): 5~40
Relative humidity: ≤80%
- Boot method: Push button to boot
- Power supply: 2×1.5VAAA alkaline batteries
- Display mode: LED
- Dimension: 58×33×33.5mm
- Weight: 35g (without batteries)
- Automatic shutdown time: 8s

QTN: 40PCS/CTN
G.W./N.W.: 4.3/1.6kg
Size: 45×26×15cm

yuwell


PORTABLE PHLEGM SUCTION UNITS

SUCTION






7E-G1



- DC power supply: AC 100-240V,50/60Hz or DC 12V
- Input power: 150VA
- Max vacuum: 85±5kPa
- Air flow rate: 27L±4/min(Device inlet)
- Operation noise: ≤65dB approximate 52db(A)
- Operating time: 30 minutes on, 30 minutes off
- Running time of battery operation: Approx. 60 minutes



7E-H1



SUCTION LINER



- Input power: 120VA
- Max vacuum: 85±5kPa
- Max air flow: 28±4L/min(Device inlet)
- Operation noise: ≤65dB(A) approximate 55db(A)
- Liquid collection jar: 1000mL, 1piece(PC material)
- Operating time: 30 minutes on, 30 minutes off

QTN: 40PCS/CTN
G.W./N.W.: 11/3kg
Size: 62.5x40x16.5cm

7E-B5

- Power supply: AC 220V, 50Hz
- Input power: 150VA
- Max vacuum: ≥80kpa
- Max air flow: 20 L/min(with jar)
- Operation noise: approx. 50dB(A) (without vacuum)
- Power: 150VA
- Liquid collection jar: 1000mL/pc
- Operating time: 30 minutes on, 30 minutes off

QTN: 4PCS/CTN
G.W./N.W.: 20.9/15.6kg
Size: 80×42×31cm





7E-A (CE/Basic sets)



Oil-free pump



Over-flow protection

- Input power: 90VA
- Max vacuum: 85±10KPA
- Max air flow: ≥20L/min (without jar)
- Operation noise: Approximately 51dB (a) (when measured at 1m from front of device)
- Liquid collection jar: 1000mL, 1piece
- Operating time: 30 minutes on, 30 minutes off

QTN: 2PCS/CTN
G.W./N.W.: 10/7.6kg
Size: 54.5×37×31cm



7E-C



Oil-free pump



Over-flow protection

- Input power: 90VA
- Max vacuum: 85±10KPA
- Max air flow: ≥20L/min (without jar)
- Operation noise: Noise: ≤65dB(A) approx.50dB(A) (without vacuum)
- Liquid collection jar: 1000mL, 1piece
- Operating time: 30 minutes on, 30 minutes off

QTN: 2PCS/CTN
G.W./N.W.: 11/8kg
Size: 52×46×34.5cm



7E-D



Over-flow protection



Battery



Oil-free pump



Vehicle use

- Power supply: AC230V, 50Hz or DC 12V
- Input power: 44VA
- Air flow rate: ≥20L/min (without jar)
- Operating time: 30 minutes on, 30 minutes off
- Running time of battery operation: Approx. 30 minutes

- Power supply of the suction apparatus and suction unit:
□AC120V±10% □AC220V±10% □AC230V±10% □50Hz □60Hz

QTN: 2PCS/CTN
G.W./N.W.: 14.7/12kg
Size: 54.5×37×31cm



7E-G



Over-flow protection



Battery



Oil-free pump



Vehicle use

- Power supply: AC230V, 50Hz or DC 12V
- Input power: 90VA
- Air flow rate: ≥20L/min (without jar)
- Operating time: 30 minutes on, 30 minutes off
- Running time of battery operation: Approx. 30 minutes

- Power supply of the suction apparatus and suction unit:
□AC120V±10% □AC220V±10% □AC230V±10% □50Hz □60Hz

QTN: 2PCS/CTN
G.W./N.W.: 14.8/12kg
Size: 52×46×34.5cm

yuwell

SUCTION APPARATUS

SUCTION






CE

NEW

7A-23A



SUCTION LINER

2L




- Input power: 180VA
- Max negative pressure: 90±10KPA
- Max air flow: 32±4L (without jar)
28±4L (with jar)
- Noise: 60dBA
- Operating time: 30 minutes on, 30 minutes off
- Liquid collection jar: 2000mL/pc, 2pieces

QTN: 1PCS/CTN
G.W./N.W.: 10.2/8.5kg
Size: 37.4×30.9×49.9cm



CE

7A-23D Electric Suction Apparatus




- Power supply: AC220V, 50Hz
- Max negative pressure: 90±10KPA
- Operation noise: ≤60 dB(A)
- Input power: 180VA
- Max air flow: ≥20L/min
- Liquid collection jar: 2500mL/pc, 2pieces

QTN: 1PCS/CTN
G.W./N.W.: 15.2/13.5kg
Size: 40.5×40×56cm



CE


7A-23B Electric Suction Apparatus



- Power supply: AC220V, 50Hz
- Max negative pressure: 90±10KPA
- Operation noise: ≤60 dB (A)
- Input power: 180VA
- Max air flow: ≥20 L/min(with jar)
40 L/min(with jar)
- Liquid collection jar: 2500mL/pc, 2pieces

QTN: 1PCS/CTN
G.W./N.W.: 19.5/17kg
Size: 44.5×40×86cm

7B Pedal Suction Apparatus



- Max negative pressure: $\geq 80\text{kPa}$
- Negative pressure back-reduction while being closed: $< 0.005\text{MPa}/10\text{min}$
- Liquid Collection Jar: 1000mL
- Pedal pressure: $\leq 245\text{N}$
- Weight: 4.2kg

QTN: 1PCS/CTN
G.W./N.W.: 4.8/3.8kg
Size: 39.5×27×31.5cm



CE


Smaller, Lighter, and Portable

No power is needed to operate this product, making it an ideal equipment for remote areas, village clinic, home care, and situations where no power is present.



CE

7B-1 Manual Suction Apparatus



- Max negative pressure: $\geq 40\text{kPa}$
- Max flow: $\geq 20\text{L}/\text{min}$
- Liquid Collection Jar: 200ml
- Weight: 202g

QTN: 5PCS/CTN
G.W.: 2.1/1.1kg
Size: 31×20×22cm



SUCTION

yuwell

ELECTRIC WHEELCHAIRS

REHABILITATION





D130AL



Fold left
and right



Electromagnetic
brake

- Aluminum alloy frame
- Fixed arm rests, fixed foot rests, foldable back
- Strong power
- Comfortable driving
- Manual and electric dual mode, easy to operate
- Outdoor type

QTN: 1PCS/CTN
G.W./N.W.: 35/29kg
Size: 108×44.5×76.5cm



Comfortable Driving



Powerful



Triple Protection



Intelligent Control



Stylish Model



D130HL



Fold left
and right



Electromagnetic
brake

- Aluminum alloy frame
- Fixed arm rests, reversible foot rests, foldable back
- Strong power
- Comfortable driving
- Manual and electric dual mode, easy to operate
- Outdoor type

QTN: 1PCS/CTN
G.W./N.W.: 31/26kg
Size: 76×41×78cm



Comfortable Driving



Powerful



Triple Protection



Intelligent Control



Stylish Model



NEW

D130TL



Fold left
and right



Electronic brake

- Aluminum alloy frame
- Fixed arm rests, elevating foot rests, foldable back
- Strong power
- Comfortable driving
- Manual and electric dual mode, easy to operate
- Outdoor type

QTN: 1PCS/CTN
G.W./N.W.: 35/29kgs
Size: 90×44×82cm



D130EL



Fold up
and down



Electromagnetic
brake




9 national
patents



Brushless
motor

- Aluminum alloy frame
- Fixed arm rests and foot rests
- Quick-release battery box
- Comfortable driving and low noise
- Manual and electric dual mode
- Easy to operate
- Outdoor type
- Brushless motor


QTN: 1PCS/CTN
G.W./N.W.: 28/22kgs
Size: 63×43×78cm




High backrest with
headrest for comfortable
riding.



Adjustable backrest and
leg rest, allowing for lying
down and resting.



Flip-up armrest, easy to
move up and down.



Compatible with backend
controller, Intelligent assi-
stance for push.



Manual and electric mode
one button switch,
easy to operate.



Extra wide tires
Anti-skid, Outstanding
gripping power.



Aerofoil foam armrests
Fatigue-reducing.



Patented supporting
backrest.
Fits the human
back curve,
extraordinary comfort.



Breathable fabrics
Lounge-style,
comfortable.



Intelligent operation
Clear display,
LED headlights,
USB port.


				
Model	D130AL	D130HL	D130TL	D130EL
Motor power	130W*2	130W*2	130W*2	130W*2
Motor type	Brush	Brush	Brush	Brushless
Lithium battery	●	●	●	●
Lead acid battery				
Electromagnetic brake	●	●		●
Parking break	●	●	●	
USB charger	●	●	●	●
Lamp				●
Quick release battery				●
Elevating mechanism			●	
Frame	Aluminum	Aluminum	Steel	Aluminum
Armrest	Fixed	Fixed	Reversible	Fixed
Footrest	Fixed	Reversible	Elevating	Reversible
Backrest	Foldable	Foldable	Foldable(90-160)	Foldable
Front wheel inch	8	8	8	8
Rear wheel inch	22	12	12	12
G.W.kg	35	31	35	28
N.W.kg	29	26	29	22
Load kg	100	100	100	100
Parking size mm	1080*445*765	760*410*780	900*440*820	650*500*770

					
	D130FL	D210BL	D210B	D210DL	D210D
	130W*2	210W*2	210W*2	210W*2	210W*2
	Brush	Brush	Brush	Brush	Brush
	●	●		●	
			●		●
	●	●	●	●	●
		●	●	●	●
	●				
	●				
	●				
	Aluminum	Aluminum	Aluminum	Steel	Steel
	Fixed	Fixed	Fixed	Reversible	Reversible
	Reversible	Reversible	Reversible	Detachable	Detachable
	Fixed	Foldable	Foldable	Foldable	Foldable
	8	8	8	8	8
	12	12	12	12	12
	27	30	40	35	45
	21	24	34	29	39
	100	100	100	100	100
	600*480*980	760*520*760	760*520*760	725*520*840	725*520*840

yuwell

MATTRESSES WITH MECHANICAL PUMP

REHABILITATION



M4 Series



Mattress coverlet with zipper



Adjustable hook



Replaceable pump filter foam



Adjustable air pump pressure


M5 Series



Mattress coverlet with zipper



Nylon + PU top cover (waterproof and moisture permeable)






Top cover with quilted cotton lining



Low pressure alarm




MATTRESSES WITH MECHANICAL PUMP				
Model	Product Image	Configuration	Selling Point	Parameter
M4-A	 QTN: 2PCS/CTN G.W./N.W: 14.2/12kg Size: 46×46×38.5cm	Pump function configuration	Low-level machinery A Two alternation Static Fireproof enclosure Stepless pressure regulation Adjustable hook design Replaceable filter cotton design Power cord 1.5 meters	Air cushion size (L x W x H): 2000x900x100mm Air pump size (L x W x H): 290x140x155mm Voltage and frequency: AC220V, 50Hz / AC110V, 50Hz Fluctuation mode: double tube alternately fluctuates Alternate time: 12min(50HZ) / 10min(60HZ) Operating pressure: 40-100mmHg Outlet flow: ≥5L/min Maximum load: 135kg Air strip specification and quantity: 4 inch pipe diameter x22 Power cord length: 1.5 meters Noise: ≤35dBA
		Bed function configuration	Head static 4 inch Inflatable tube Nylon/PVC gray Upper cover-polyester red Base cushion Nylon/PVC dark gray Base cushion fixing strap All—around zipper	
M5-A	 QTN: 2PCS/CTN G.W./N.W: 18/15.2kg Size: 53×49×40cm	Pump function configuration	Two alternation Static Fireproof enclosure Stepless pressure regulation Adjustable hook design Replaceable filter cotton design Power cord 4.5 meters Handle Low—pressure warning light	Air cushion size (L x W x H): 2000x900x125mm Air pump size (L x W x H): 300x140x190mm Voltage and frequency: AC220V, 50Hz / AC110V, 60Hz Fluctuation mode: double tube alternately fluctuates Alternate time: 12min(50HZ) / 10min(60HZ) Operating pressure: 30-80mmHg Outlet flow: ≥5L/min Maximum load: 160kg Air strip specification and quantity: 5 inch pipe diameter x18 Power cord length: 4.5 meters Noise: ≤30dBA
		Bed function configuration	Head static 5 inch Inflatable tube Nylon/PVC gray Upper cover-polyester + car cotton red Base cushion Nylon/PVC dark gray Base cushion fixing strap All—around zipper	
M5-B	 QTN: 2PCS/CTN G.W./N.W: 18/15.2kg Size: 53×49×40cm	Pump function configuration	Two alternation Static Fireproof enclosure Stepless pressure regulation Adjustable hook design Replaceable filter cotton design Power cord 4.5 meters Handle Low—pressure warning light	Air cushion size (L x W x H): 2000x900x125mm Air pump size (L x W x H): 300x140x190mm Voltage and frequency: AC230V, 50Hz / AC110V, 60Hz Fluctuation mode: double tube alternately fluctuates Alternate time: 12min(50HZ) / 10min(60HZ) Operating pressure: 30-80mmHg Outlet flow: ≥5L/min Maximum load: 160kg Air strip specification and quantity: 5 inch pipe diameter x18 Power cord length: 4.5 meters Noise: ≤30dBA
		Bed function configuration	Head static Micro-jet 5 inch Inflatable tube Nylon/PVC gray Upper cover-polyester + car cotton red Base cushion Nylon/PVC dark gray Base cushion fixing strap All—around zipper	
M5-C	 QTN: 2PCS/CTN G.W./N.W: 18/15.2kg Size: 53×49×40cm	Pump function configuration	Two alternation Static Fireproof enclosure Stepless pressure regulation Adjustable hook design Replaceable filter cotton design Power cord 4.5 meters Handle Low—pressure warning light	Air cushion size (L x W x H): 2000x900x125mm Air pump size (L x W x H): 300x140x190mm Voltage and frequency: AC220V, 50Hz / AC110V, 60Hz Fluctuation mode: double tube alternately fluctuates Alternate time: 12min(50HZ) / 10min(60HZ) Operating pressure: 30-80mmHg Outlet flow: ≥5L/min Maximum load: 160kg Air strip specification and quantity: 5 inch pipe diameter x18 Power cord length: 4.5 meters Noise: ≤30dBA
		Bed function configuration	Head static CPR transport 5 inch Inflatable tube Nylon/PVC gray Upper cover-polyester + car cotton red Base cushion Nylon/PVC dark gray Base cushion fixing strap All—around zipper	

M8-A



Static mode
Convenient for patients to turn over and nursing.



Alternate mode
Alternately simulate natural sleep rhythm, promote blood circulation and avoid local ischemia caused by long-term compression of tissues.



Head restraint function
The patient's head is located in the non fluctuating front four tube headrest airbag, which makes the patient sleep more comfortably.



CPR quick release valve
When cardiopulmonary resuscitation is needed, be sure to open CPR to quickly deflate the mattress.



Braided band
The shuttleless Loom makes braided band, which is more durable.



Quick replacement
It is convenient for users to quickly remove and replace the airbag.



Two way zipper
It can be assembled and disassemble at will, which is very convenient for cleaning when mixed disinfection is carried out in medical institutions.



Mattress micro jet
Reduce the incidence of pressure ulcer, avoid the accumulation of moisture and heat in the bed, and achieve the effect of dryness and comfort.




Nylon Pu Bedsread
Feel comfortable and silky. Easy to clean and avoid dirt sticking.



- Adjustable weight
- Low noise
- Dual mode
- Alarm Mute Button

Panel information:

-  Knob adjustment pressure
-  Dual mode
-  Pressure indicator




MATTRESSES WITH MECHANICAL PUMP				
Model	Product Image	Configuration	Multi-function	Parameter
M8-A	 QTN: 1PC/CTN GW./NW.: 12.6kg/11.4kg Size: 46.5cm×37cm×43cm	Pump function configuration	Alternation mode Static mode Fireproof enclosure Stepless pressure regulation Replaceable air filter cotton 4.5 meters power cord Pump handle Low-pressure alarm Alarm Mute	Air cushion size(L x W x H): 2000x900x200mm Air pump size (L x W x H) 290x140x190mm Voltage and frequency: AC220V~240V,50Hz/AC110V~120V,60Hz Fluctuation mode: double tube alternately fluctuates Alternate time: 12min(50Hz)/10min(60Hz) Operation pressure: 30-80mmHg Outlet flow: 8L/min Maximum load: 220kg Air strip specification and quantity: 8 inch pipe diameter × 20 Power cord length: 4.5meters Noise: ≤ 30dBA
		Bed function configuration	Static headrest Micro-jet airbag 8 inch mattress Inflatable tube Nylon/PVC blue Upper cover Nylon/PU blue Base cushion Nylon/PVC blue Base cushion fixing strap All-around zipper	



REHABILITATION

YUWELL IOT PLATFORM



GPRS / WIFI Bluetooth®



BreathCare PAP

GPRS



Oxygen Concentrator

Bluetooth®



Blood Pressure Monitor

Bluetooth®



Blood Glucose Meter

Bluetooth®



Thermometer

Bluetooth®



Oximeter